

Moi-Stik™

Oral hygiene+

MOI-STIK™

Offering an easy to use, gentle but effective way to help manage oral hygiene and Dry Mouth on the ward and in the care home environment.

The Moi-Stik is supplied in handy sealed foil packs with three soft, moist, cotton tipped applicators.



The Moi-Stick are supplied in a box of 50 foil packs.

The large, soft cotton applicators (the tips of which are both adhered and tightly wound around the 4" shaft) are irradiated before being moistened with Saliveze. This pH neutral, alcohol and sugar free aqueous solution contains important electrolytes, and is one of the most widely prescribed Dry Mouth solutions in the community.

Perfect for any ward or care environment, simply move the Moi-Stik around all intra-oral surfaces, floor of the mouth, tongue and lips. This helps to gently cleanse the mouth, whilst offering important lubrication via the Saliveze solution.

Moi-Stik can help a patient's QOL through the relief of Dry Mouth symptoms by:

- Lubricating the oral mucosa
- Maintaining a neutral pH
- Aiding speaking and swallowing
- Replacing salivary electrolytes
- Helping reduce tooth decay
- Aiding restful sleep

Saliveze Dry Mouth Spray
50ml pump spray - On prescription



Ordering information

Code	Description	NHS Supply Chain
MS-50	Moi-Stik - Box of 50 sachets	NA - direct only
MS-50-4	Moi-Stik - Case of 4 x Boxes of 50 sachets.	IEZ259
MS-50-10	Moi-Stik - Case of 10 x Boxes of 50 sachets.	IEZ255

For more information or a sample please contact us
T: 01531 631105 E: info@wyvernmedical.co.uk
W: www.wyvernmedical.co.uk